

Medical Device

Collagen			
Intracranial Hypertension	Size	Applied Fields	Note
CollaDental Matrix (Purified Porcine type I Collagen)	8 x 15mm 8 x 25mm 15 x 25mm	Intended for using in dental surgical procedures as a resorbable material for placement in the area of dental implants, bone defects, extraction socket or ridge reconstruction.	
CollaDental Graft (60% HAP(Hydroxyapatite), 40% β -TCP(Tricalcium Phosphate) and Purified Porcine type I Collagen)	0.25g 0.5g	<ul style="list-style-type: none"> ✧ Augmentation or reconstructive treatment of the alveolar ridge. ✧ Filling of periodontal defects. ✧ Filling of bone defects after root resection, apicoectomy and cystectomy. ✧ Filling of extraction sockets to enhance preservation of the alveolar ridge. ✧ Elevation of the maxillary sinus floor. ✧ Filling of periodontal defects in conjunction with products intended for Guided Tissue. Regeneration(GTR) and Guided Bone Regeneration(GBR). ✧ Filling of peri-implant defects in conjunction with products intended for Guided Bone Regeneration(GBR). 	
CollaDental Barrier (Purified Porcine type I Collagen)	15 x 20mm 20 x 30mm 30 x 40mm	<ul style="list-style-type: none"> ✧ Localized ridge reconstruction for later implantation. ✧ Alveolar ridge reconstruction for prosthetic treatment. ✧ Filling of bone defects after root resection, cystectomy, removal of retained teeth. ✧ Guided bone refeneration in dehiscence defects ✧ Guided tissue regeneration procedures in periodontal defects. 	