

Medical Device

Sinbone			
Intracranial Hypertension	Size	Applied Fields	Note
Sinbone (Ceramic Bone Graft) 60% HAP(Hydroxyapatite) from bovine bone(Immunogenic-free) and 40% β -TCP(Tricalcium Phosphate)	5x5x5mm (8's Cubic/1cm ³ /Box)	Orthopaedic Surgery, Neurology Surgery, Plastic Surgery and Otolaryngology. Indicated for: ✧ Bone defects after removal of autologous bone grafts. ✧ Bone defects after removal of tumors. Pseudarthrosis. ✧ Spondylodesis (with only low axial compression of the implant). ✧ Acetabular reconstruction in child . ✧ Bone defects after revision of loosened prosthetic component.	ISO 9001 ISO 13485
	5x5x20mm (4's Strip/2cm ³ /Box)		
	5x5x20mm (10's Strip/5cm ³ /Box)		
	0.25mm-0.5mm (0.25g Granule/Box)		
	0.25mm-0.5mm (0.5g Granule/Box)		
	1.0mm-2.0mm (0.5g Granule/Box)		
	1.0mm-2.0mm (1.5g Granule/Box)		
	0.5mm-1.0mm (1cm ³ Granule/Box)		
	0.5mm-1.0mm (2cm ³ Granule/Box)		
	0.5mm-1.0mm (7cm ³ Granule/Box)		